

Autoclavable Biohazard Bag

PMSSteripack[™]

PMSSteripack autoclavable biohazard bags are ideal for processing of biohazardous materials which require autoclaving prior to disposal. They provide a safe handling and efficient solution for your biohazard waste management.

High integrity polypropylene bag withstand 141°C (285°F) maximum autoclave temperature and is designed to resist puncture, tears and leaks.

PMSSteripack autoclavable biohazard bags are available in different sizes and all bags are imprinted with the universal biohazard symbol.



Features

- Red, strong autoclavable bag for biohazard waste
- Designed for STEAM sterilization at max. 141°C
- Puncture, tear and leak resistance
- Will not soften or melt during sterilization
- Clearly imprinted biohazard symbol

Ref.	Dimension (L x H)	Unit
OT5370	530mm x 700mm	200pcs/box
OT6380	630mm x 800mm	120pcs/box
OT6389	630mm x 890mm	120pcs/box
OT80110	800mm x 1100mm	100pcs/box
OT87110	870mm x 1100mm	100pcs/box
	OT5370 OT6380 OT6389 OT80110	OT5370 530mm x 700mm OT6380 630mm x 800mm OT6389 630mm x 890mm OT80110 800mm x 1100mm

