

# TECHNICAL DATA SHEET

## AUTOCLAVE TAPE (UNPRINTED)

<b>PRODUCT CODE</b>	:	<b>AT1950 – AT2550</b>
<b>RELATED STANDARDS</b>	:	CE marking is not applicable for autoclave tape products, because this product is not considered as medical device under the clause of the Medical Device Regulation (Eu) 2017/745
<b>INTENDED USE</b>	:	Unprinted Autoclave Tapes are designed to close sterilization packs.
<b>OPERATION CONDITIONS</b>	:	Unprinted autoclave tape are used in all types of sterilizers except for Plasma ( VH2O2 ) Type.
<b>CHEMICAL COMPOSITION</b>	:	It doesn't contain any chemical component except Hot Melt Adhesive. The product is not classified as dangerous.
<b>PRODUCT SPECIFICATIONS</b>	:	

AUTOCLAVE TAPE	SPECIFICATION	UNIT	STANDARD	TYPICAL VALUE
WIDTH		mm	-	19±1
				25±1
LENGTH		m	-	50± % 1
THICKNESS		µm	-	100 ±5
ADHESIVE		g/m <sup>2</sup>	-	35 ±5
INSTANT TACK		N/19mm	ASTM D 3121	Min.1.0
PEEL ADHESION		N/19mm	ASTM D 3121	Min.1.3

**PACKAGING** : Inner Packaging : Each product is wrapped with polyester packaging.

Package Contents: 19 mm x 50 m, 60 pieces in carton box  
Package Contents: 25 mm x 50 m, 36 pieces in carton box



PMSSteriT<sup>™</sup>

PMS Tıbbi Cihazlar Teknolojisi San. ve Tic. A.Ş.

pms@pmsmedikal.com