

TECHNICAL DATA SHEET

STERILIZATION BASKETS

- PRODUCT CODE** : **TM**
- INTENDED USE** : Crimp Wire Baskets are ideal for hospitals and clinics to hold several instruments at once during the washing, disinfection, sterilization and storage.
- OPERATION CONDITIONS** : Crimp Wire Baskets are used for steam, ethylene oxide, formaldehyde, gas plazma (VH2O2) sterlization.The sterilization conditons should be determined by the end user regarding to material to be sterilized.
- PRODUCT SPECIFICATIONS** : PMS crimp wire baskets are made of high quality AISI 304 stainless steel and are available in different sizes and designs.

PRODUCT SIZE & PACKAGING:

Code	Description	Size	Unit/Box
TM.100.030	Crimp Wire Basket 1/1	540x255x30 mm	1 pcs
TM.100.050	Crimp Wire Basket 1/1	540x255x50 mm	1 pcs
TM.100.070	Crimp Wire Basket 1/1	540x255x70 mm	1 pcs
TM.100.100	Crimp Wire Basket 1/1	540x255x100 mm	1 pcs
TM.105.030	Crimp Wire Basket 1/1	485x255x30 mm	1 pcs
TM.105.050	Crimp Wire Basket 1/1	485x255x50 mm	1 pcs
TM.105.070	Crimp Wire Basket 1/1	485x255x70 mm	1 pcs
TM.105.100	Crimp Wire Basket 1/1	485x255x100 mm	1 pcs
TM.140.030	Crimp Wire Basket 3/4	405x255x30 mm	1 pcs
TM.140.050	Crimp Wire Basket 3/4	405x255x50 mm	1 pcs
TM.140.070	Crimp Wire Basket 3/4	405x255x70 mm	1 pcs
TM.140.100	Crimp Wire Basket 3/4	405x255x100 mm	1 pcs
TM.160.030	Crimp Wire Basket 1/2	255x245x30 mm	1 pcs
TM.160.050	Crimp Wire Basket 1/2	255x245x50 mm	1 pcs
TM.160.070	Crimp Wire Basket 1/2	255x245x70 mm	1 pcs
TM.160.100	Crimp Wire Basket 1/2	255x245x100 mm	1 pcs
TM.101.030	Lid for Baskets	510x240 mm	1 pcs
TM.106.030	Lid for Baskets	455x240 mm	1 pcs
TM.141.030	Lid for Baskets	375x240 mm	1 pcs
TM.161.030	Lid for Baskets	240x215 mm	1 pcs
TM.420.010	Basket Tablet	240x240x32 mm	1 pcs
TM.420.020	Basket Tablet	425x240x32 mm	1 pcs



PMS Steripack®