## **TECHNICAL DATA SHEET**

## **AUTOCLAVABLE BIOHAZARD BAG**

PRODUCT CODE : OT

**RELATED** : Autoclavable Biohazard Bag meets the requirements of the

STANDARDS Council Directive 94/62/EC on Packaging.

INTENDED USE : Autoclavable Biohazard Bags are used for the sterilization of

medical wastes before disposal.

OPERATION : Autoclavable Biohazard Bags are designed for use in steam sterilization methods. The recommended operation temperature

is max. 141 °C. Other sterilization conditions should be determined by the end user regarding to material to be sterilized. During sterilization, bag should be left open for an

effective sterilization.

## **PRODUCT SPECIFICATIONS**

PROPERTIES	METHOD	UNIT	VALUE	
THICKNESS (one side)	ISO 534	μm	50	
TENSILE STRENGTH	ISO 527-3	N/15 mm	MD	14
			CD	12,5
ELONGATION AT BREAK	ISO 527-3	%	MD	320
			CD	120
TEAR FORCE	ASTM D 1922	mN	3300	

**PACKAGING** : Autoclavable Biohazard Bags are packaged as the listed below.

Product Code	<b>Dimensions</b>	<u>Unit</u>
OT5370	53 cm x 70 cm	300 pcs/box
OT6389	63 cm x 89 cm	200 pcs/box
OT8790	87 cm x 90 cm	200 pcs/box
OT87110	87 cm x 110 cm	120 pcs/box
OT80110	80 cm x 110 cm	100 pcs/box

